

产品说明书

品名: HA Peptide

货号: HA-P02 / HA-P05 / HA-P25 / HA-P100

规格: 2/5/25/100mg

Alias: HA Tag; HA Synthetic Peptide; HA tag peptide;

➤Description:

A fusion tag, called HA-tag and consisting of 9 amino acids (YPYDVPDYA) , was specifically designed for immunoaffinity chromatography. It allows elution under non-denaturing conditions.

➤Sequence (one-letter code): YPYDVPDYA

➤Sequence (three-letter code): {Tyr}{Pro}{Tyr}{ASP}{Val}{Pro}{ASP}{Tyr}{Ala}

➤Formula: C₅₃H₆₇N₉O₁₇

➤Molecular Weight: 1,102

➤Purity: > 95% by HPLC

➤Storage: Store at -20°C for at least 12 months

For Laboratory Research Use Only